





Maintain the pushed locator, insert the contact, release then place the bare thread.

A pressure on sleeves and the contact is set according to our crimping specifications.

## **MULTI CONTACT HAND TOOL**

- Our PCR Crimp Tools are equipped with locators which optimize contact installation and wire insertion. Therefore they simplify the achievement of crimp cycle.
- Their handles reduce crimping effort. The symmetric crimp dies and locator can be used by right and left-handed users. These tools are compact and light. They allow an accurate crimping for the professional users in search of constant reliability and quality.
- Lifetime tests exceeding 50 000 cycles.



<b>FAMILY</b>	TERMINALS	<b>CRIMP TOOL</b>	WIRE SECTION (mm <sup>2</sup> )	AWG
PHD	SPHD-001T-P0.5	PCR-04AB	0.13 - 0.33	26 – 22
	SPHD-002T-P0.5		0.08 - 0.21	28 - 24
PA	SPA-001T-P0.5		0.13 - 0.33	26 – 22
PAL	SPAL-001T-P0.5		0.13 - 0.33	26 – 22
	SPAL-002T-P0.5		0.08 - 0.21	28 – 24



## Features:

Length: 190 mm
Height: 74 mm
Width: 21 mm

Weight:

280 g (hand tool only)

302 g (hand tool + die)

317 g (hand tool + die + locator)

